



# Accuracy of a Single, Heparin-Calibrated Anti-Xa Assay for the Measurement of Rivaroxaban, Apixaban, and Edoxaban Drug Concentrations: A Prospective Cross-Sectional Study

Tamana Meihandoest<sup>1,2</sup>, Jan-Dirk Studt<sup>3,4</sup>, Adriana Mendez<sup>5</sup>, Lorenzo Alberio<sup>6</sup>, Pierre Fontana<sup>7</sup>, Walter A. Wuillemin<sup>8</sup>, Adrian Schmidt<sup>9</sup>, Lukas Graf<sup>10</sup>, Bernhard Gerber<sup>11,4</sup>, Ursula Amstutz<sup>2</sup>, Cedric Bovet<sup>2</sup>, Thomas C. Sauter<sup>12</sup>, Lars M. Asmis<sup>13</sup> and Michael Nagler<sup>2\*</sup>

<sup>1</sup> Department of Epidemiology, Maastricht University, Maastricht, Netherlands, <sup>2</sup> Department of Clinical Chemistry, Inselspital, Bern University Hospital and University of Bern, Bern, Switzerland, <sup>3</sup> Department of Medical Oncology and Hematology, University Hospital Zurich, Zurich, Switzerland, <sup>4</sup> University of Zurich, Zurich, Switzerland, <sup>5</sup> Institute for Laboratory Medicine, Kantonsspital Aarau AG, Aarau, Switzerland, <sup>6</sup> Service and Central Laboratory of Hematology, CHUV, Lausanne University Hospital, Lausanne, Switzerland, <sup>7</sup> Division of Angiology and Hemostasis, Geneva University Hospitals, Geneva, Switzerland, <sup>8</sup> Division of Hematology and Central Hematology Laboratory, Cantonal Hospital of Lucerne and University of Bern, Bern, Switzerland, <sup>9</sup> Institute of Laboratory Medicine and Clinic of Medical Oncology and Hematology, City Hospital Zurich, Zurich, Switzerland, <sup>10</sup> Center for Laboratory Medicine, St. Gallen, Switzerland, <sup>11</sup> Clinic of Hematology, Oncology Institute of Southern Switzerland, Bellinzona, Switzerland, <sup>12</sup> Department of Emergency Medicine, Inselspital, Bern University Hospital, Bern, Switzerland, <sup>13</sup> Centre for Perioperative Thrombosis and Haemostasis, Zurich, Switzerland

**Background:** Applying a single anti-Xa assay, calibrated to unfractionated heparin to measure rivaroxaban, apixaban, and edoxaban would simplify laboratory procedures and save healthcare costs.

**Aim:** We hypothesized that a heparin-calibrated anti-Xa assay would accurately measure rivaroxaban, apixaban, and edoxaban drug concentrations and correctly predict clinically relevant drug levels.

**Methods:** This analysis is part of the Simple-Xa study, a prospective multicenter cross-sectional study conducted in clinical practice. Patients treated with rivaroxaban, apixaban, or edoxaban were included. Anti-Xa activity was measured using the Siemens INNOVANCE<sup>®</sup> Heparin assay. Drug concentrations were determined using ultra-high performance liquid chromatography-tandem mass spectrometry (LC-MS/MS). Cut-off levels were determined in a derivation dataset (50% of patients) and sensitivities and specificities were calculated in a verification dataset (50% of patients).

**Results:** Overall, 845 patients were available for analysis. Correlation coefficients ( $r_s$ ) between the heparin-calibrated anti-Xa assay and drug concentrations were 0.97 (95% Cl 0.97, 0.98) for rivaroxaban, 0.96 (0.96, 0.97) for apixaban, and 0.96 (0.94, 0.99) for edoxaban. The area under the receiver operating characteristics curve (ROC) was 0.99

#### OPEN ACCESS

#### Edited by:

Roberto Pola, Università Cattolica del Sacro Cuore, Italy

#### Reviewed by:

Francois Mullier, CHU Dinant Godinne UCL Namur, Belgium Heiko Rühl, University Hospital Bonn, Germany

#### \*Correspondence:

Michael Nagler Michael.Nagler@insel.ch

#### Specialty section:

This article was submitted to Thrombosis, a section of the journal Frontiers in Cardiovascular Medicine

> Received: 18 November 2021 Accepted: 22 February 2022 Published: 17 March 2022

#### Citation:

Meihandoest T, Studt J-D, Mendez A, Alberio L, Fontana P, Wuillemin WA, Schmidt A, Graf L, Gerber B, Amstutz U, Bovet C, Sauter TC, Asmis LM and Nagler M (2022) Accuracy of a Single, Heparin-Calibrated Anti-Xa Assay for the Measurement of Rivaroxaban, Apixaban, and Edoxaban Drug Concentrations: A Prospective Cross-Sectional Study. Front. Cardiovasc. Med. 9:817826. doi: 10.3389/fcvm.2022.817826

1

for all clinically relevant drug concentrations. In the verification dataset, the sensitivity was 94.2% (95% Cl 90.8–96.6) for 30  $\mu$ g L<sup>-1</sup>, 95.8% (92.4–98.0) for 50  $\mu$ g L<sup>-1</sup>, and 98.7% (95.5–99.9) for 100  $\mu$ g L<sup>-1</sup>. Specificities were 86.3% (79.2–91.7), 89.8% (84.5–93.7), and 88.7% (84.2–92.2), respectively.

**Conclusion:** In a large prospective study in clinical practice, a strong correlation of heparin-calibrated anti-Xa measurements with LC-MS/MS results was observed and clinically relevant drug concentrations were predicted correctly.

Keywords: diagnostic accuracy, anti-Xa assay, laboratory monitoring, direct oral anticoagulants, rivaroxaban

# HIGHLIGHTS

#### What is known about this topic?

- Applying a single anti-Xa assay to measure rivaroxaban, apixaban, and edoxaban would simplify laboratory procedures and save healthcare costs.
- It remains unclear if this can be achieved using a single anti-Xa assay, calibrated to unfractionated heparin.

#### What does this paper add?

- We conducted a prospective multicenter cross-sectional study including 845 patients taking rivaroxaban, apixaban, or edoxaban in clinical practice.
- The association between heparin-calibrated anti-Xa measurements and LC-MS/MS results was strong for all drugs.
- Clinically relevant drug levels were predicted correctly.

# INTRODUCTION

The proportion of patients taking direct oral anticoagulants (DOAC) for the prevention and treatment of thromboembolic diseases is rapidly increasing (1, 2). These patients occasionally face clinical situations with a high bleeding risk such as accidents, urgent surgery, and thrombolysis because of acute stroke (3-5). Besides, relevant DOAC drug concentrations contribute to massive bleeding of any cause (3, 6, 7). In addition, unresponsive or demented patients present to the emergency department, and knowledge of anticoagulant treatment is essential in the management. Rapid determination of DOAC drug levels in these situations supports clinical decisions regarding reversal agents, and deferral of interventions (8-10). Additionally, accumulation in the case of renal and/or hepatic failure or even overdosing can be detected (11, 12). Thus, determination of DOAC drug levels in special clinical situations is recommended by major scientific societies such as the International Society on Thrombosis and Haemostasis (9). Furthermore, it might potentially save health care costs associated with the clinical situations mentioned above (13).

Ideally, a simple laboratory test that accurately determines DOAC plasma levels would be available and implemented in various healthcare settings in a 24/7 service (3). Routine coagulation tests are neither sensitive nor specific in the detection of DOAC (10, 14–16). Various anti-Xa assays are available that measure rivaroxaban, apixaban, or edoxaban concentrations using drug-specific calibration curves (17). However, these tests are still not widely implemented because it is laborious and expensive to provide three different tests (3). We and other authors hypothesized that a single heparincalibrated anti-Xa assay would be sufficient to accurately and efficiently determine rivaroxaban, apixaban, and edoxaban drug levels (18-20). An essential advantage of this method is that heparin-calibrated assays are already available even in smaller laboratories. Thus, implementing a single-calibration anti-Xa assay for unfractionated heparin, low molecular weight heparin, rivaroxaban, apixaban, and edoxaban would improve laboratory procedures but also the care of patients treated with DOAC. Recently, we demonstrated that a universal, LMWH-calibrated assay can accurately measure DOAC drug levels (21). However, it remains unclear if this can be achieved using a universal heparincalibrator intended to detect unfractionated heparin and low molecular weight heparins.

With the present multicenter cross-sectional study, we aimed to assess whether the Siemens INNOVANCE<sup>®</sup> Heparin anti-Xa assay would accurately measure rivaroxaban, apixaban, and edoxaban drug concentrations and correctly predict clinically relevant drug levels. The primary focus was on the clinically significant concentration range between 0 and 300  $\mu$ g L<sup>-1</sup>.

# MATERIALS AND METHODS

# Study Design, Setting, and Population

We conducted a prospective multicenter cross-sectional study including patients in nine hemostasis laboratories affiliated to Swiss tertiary hospitals. Patients treated with rivaroxaban, apixaban, or edoxaban in clinical practice were included between 2018 and 2019 (**Figure 1**, CONSORT flow diagram). Inclusion criteria were (a) 18 years or older, (b) use of rivaroxaban, apixaban, or edoxaban (c) DOAC drug-level requested, and (d) signed general informed consent, if required by local authorities. Exclusion criteria were (a) refused general informed consent, (b) use of heparin, (c) preanalytical issues, (d) use of more than one DOAC, and (e) insufficient sample material. Samples were collected regardless of the time of last drug intake, covering the full range of drug levels observed in clinical practice. Ultra-high performance liquid chromatography-tandem mass spectrometry (LC-MS/MS) measuring rivaroxaban, apixaban, and edoxaban



was used as a reference standard (16, 21). The study was approved by the local ethics committees and all hospitals gave local feasibility approval. If required, patients signed a general consent to use their samples and data before enrolment at the respective study center. The study was conducted in accordance with the declaration of Helsinki.

# **Data Collection and Handling of Samples**

Patient characteristics including age, sex, and the DOAC used (rivaroxaban, apixaban, or edoxaban) were collected in a secured REDCap database. Protocols were implemented at all institutions detailing blood-drawing procedures to ensure adequate preanalytical conditions (22). Venous blood samples were drawn in plastic syringes containing 1 mL trisodium citrate (0.106 mol  $L^{-1}$ ) per mL of blood (i.e., S-Monovette<sup>®</sup> Citrate, Sarstedt, Nümbrecht, Germany). Samples were centrifuged according to an established scheme (1,500 g for 10–15 min, or 3,137 g for 7 min), and aliquots were frozen immediately and stored at  $-80^{\circ}$ C until transportation (22). Samples were shipped on dry ice in one batch per site to the central laboratory and delivered within 3 to 4 h. Samples were kept frozen until the determination of laboratory tests without any freeze-thaw cycle. Laboratory test results were exported automatically to avoid typing errors.

# **Determination of the Anti-Xa Activity**

The Siemens INNOVANCE® Heparin anti-Xa assay was selected for determining the anti-Xa activity. This one-stage chromogenic assay is designed to measure the drug level of unfractionated heparin (UFH) as well as low molecular weight heparin (LMWH) (INNOVANCE® Heparin anti-Xa assay application sheet; Siemens Healthineers, Marburg, Germany). The reagent contains dextran sulfate but not exogenous antithrombin. The instructions provided by the manufacturer were strictly followed. In brief, a five-step calibration curve was applied using the calibrators provided by the manufacturer (0.0, 0.42, 0.86, 1.26, and 1.67 U/ml). Samples were rapidly thawed and gently mixed at 37°C. Patient plasma was pre-diluted (1:2) and added to the reagent containing coagulation factor Xa and a chromogenic substrate. The formation of paranitroaniline was quantified



optically at a wavelength of 405 nm. The measurements were conducted in batches on an Atellica COAG 360 analyzer (Siemens Healthineers, Marburg, Germany) (23). Measurements were performed blinded to the LC-MS/MS test results.

# Determination of the Drug Concentration by Liquid Chromatography-Tandem Mass Spectrometry

Liquid chromatography-tandem mass spectrometry was used for the quantification of rivaroxaban, apixaban, edoxaban, and M4 metabolite of edoxaban. For protein precipitation and analyte extraction, plasma acetonitrile:water 1:1 (v/v), extraction buffer (MassTox TDM Series A, Chromsystems, Gräfelfing, Germany), and precipitation reagent (MassTox TDM Series, Chromsystems, Gräfelfing, Germany) containing the isotope labeled internal standards (13C6 rivaroxaban, 13CD3 apixaban, 13CD2 edoxaban; provided by the manufacturers) were added to the plasma. Afterward, the samples were vortexed and centrifuged at 14,000 rcf and 20°C for 4 min. The supernatant was diluted with water:methanol 8:2 (v/v) and stored at 10°C until analysis. The calibrators and QCs were prepared in pooled plasma (Innovative Research, Novi, MI, United States). The extracted samples were analyzed using reversed-phase chromatography (Cortecs UPLC C18 column,  $2.1 \times 75$  mm,  $1.7 \,\mu$ m, Waters) on a triple quadrupole mass spectrometer (Xevo TQ-S, Waters, Milford, CT, United States) coupled to a UPLC Acquity I-Class system (Waters, Milford, CT, United States). Edoxaban M4 concentration was summed up with the edoxaban for further analysis.

# **Statistical Analysis**

Variables were described using proportions and percentages or median and interquartile range (IQR) as appropriate. The normality of the data was assessed visually and using a Q-Q plot. The accuracy of the anti-Xa assay was determined by calculating the Spearman's correlation coefficient in relation to the plasma concentration as measured by LC-MS/MS. A correlation coefficient of  $rs \ge 0.95$  was considered as accurate (alternative hypothesis), and a correlation of  $rs \le 0.6$  was regarded as inadequate (null hypothesis). The Deming regression was used to describe the linear relationship and a modified Bland-Altman plot (ratios were used due to different scales) was created

TABLE 1   Characteristics of patients treated with rivaroxaban,	apixaban or
edoxaban ( $n = 845$ ).	

	Patients treated with					
	Rivaroxaban	Apixaban	Edoxaban	All	Missing data	
Patients (n/%)	374 (44.3)	405 (47.9)	66 (7.8)	845 (100)	0	
<b>Age</b> (years; median/IQR)	74 (63-83)	78 (63-83)	75 (58-82)	76 (66-82)	96	
Sex (n/%)					8	
Male	209 (56.0)	233 (58.4)	36 (55.4)	478 (57.1)	-	
Female	164 (44.0)	166 (41.6)	29 (44.6)	359 (42.9)	-	

N, number. IQR, interquartile range.

to observe a potential bias over the spectrum of measurements (9). Systematic differences were analyzed by calculating the mean difference and the SD to compute 95% limits of agreement for every level of measurements (average difference  $\pm$  1.96 standard deviation of the difference) (24). To assess the diagnostic accuracy of the anti-Xa assay, we determined the sensitivity and specificity of detecting 30, 50, and 100  $\mu$ g/L, representing clinically relevant drug levels. The dataset was randomly split in half, and the cutoffs of the new test were obtained by a ROC curve analysis in the derivation dataset. As an internal validation, we repeated

the ROC curve analysis in the verification dataset and calculated sensitivities and specificities regarding clinically relevant drug levels. An area under the receiver-operating characteristics (ROC) curve  $\geq 0.95$  and sensitivities/specificities of at least 90% were regarded as adequate. A power analysis for a one-sample correlation test was conducted with a power of 0.9 and an alpha of 0.05. Since many patients were taking rivaroxaban and apixaban, and only a few edoxaban, 932 patients were included until data saturation was reached. Sensitivity analyses considering samples below  $3.34 \text{ U ml}^{-1}$  only were conducted. All statistical analyses were performed using RStudio (1.3. 1093-1); figures were created using Prism 8 (GraphPad Software, Inc., La Jolla, CA, United States).

## RESULTS

### **Patient Characteristics**

Overall, 932 patients were included in this prospective multicenter cross-sectional study; a detailed flow chart is given in **Figure 1**. From this study population, 35 patients were excluded because of heparin use, five patients due to preanalytical issues, two patients using more than one DOAC, and 45 patients due to insufficient sample material. Eventually, samples of 845 patients were used for the current analysis. Of these, 374 patients used rivaroxaban, 405 apixaban, and 66 edoxaban. The median



practice. Ultra-performance liquid chromatography-tandem mass spectrometry (LC-MS/MS) was used to determine drug levels. Spearman's correlation coefficient (r<sub>s</sub>) was 0.97 for rivaroxaban (95% CI 0.97–0.98), 0.96 for apixaban (0.95–0.97), and 0.96 for edoxaban (0.94–0.99). The overall r<sub>s</sub> was 0.97 (95% CI, 0.97–0.98).

TABLE 2 | Accuracy of universal, heparin-calibrated anti-Xa assay with regard to drug concentration in 845 patients taking rivaroxaban, apixaban, or edoxaban in clinical practice.

	Overall	Rivaroxaban	Apixaban	Edoxaban		
Spearman's correlation coefficient (95% Cl)	0.97 (0.97, 0.98)	0.97 (0.97, 0.98)	0.96 (0.96, 0.97)	0.96 (0.94, 0.99)		
Deming regression Slope (95% Cl)	0.008 (0.006, 0.009)	0.008 (0.007, 0.009)	0.008 (0.007, 0.009)	0.006 (0.005, 0.007)		
Y-intercept (95% CI)	0.29 (0.15, 0.42)	0.29 (0.16, 0.40)	0.14 (0.05, 0.23)	0.10 (0.015, 0.18)		

Ultra-performance liquid chromatography-tandem mass spectrometry (LC-MS/MS) was used to determine drug levels. The Spearman's correlation coefficient is given ( $r_s$ ), and the coefficients of the Deming regression.



age was 76 years old (IQR, 66 to 82 years), and 42.9% of the patients were female, see Table 1.

# Association Between Anti-Xa Activity and Drug Concentration

The association between anti-Xa activity and drug concentrations as measured by LC-MS/MS is shown in **Figure 2**, and results of correlation analyses are provided in **Table 2**. The overall Spearman's correlation coefficient ( $r_s$ ) was 0.97 (95% confidence interval, CI, 0.97 to 0.98). Regarding the individual drugs,  $r_s$ was 0.97 (95% CI 0.97–0.98) for rivaroxaban, 0.96 (0.95–0.97) for apixaban, and 0.96 (0.94–0.99) for edoxaban. The overall slope of the regression equation was 0.008 (95% CI 0.007– 0.009), 0.008 (0.007–0.009) for rivaroxaban, 0.009 (0.008–0.010) for apixaban, and 0.006 (0.005–0.008) for edoxaban. The overall Y-intercept was 0.21 (95% CI 0.13–0.28), 0.29 (0.16–0.40) for rivaroxaban, 0.12 (0.016–0.022) for apixaban, and 0.10 (0.007–0.16) for edoxaban. The association between anti-Xa activity and drug concentrations over the range of measurements is shown in a modified Bland-Altman plot, see **Supplementary Figure 1**. The overall bias was 0.01, with a lower limit of agreement of -0.01, and an upper limit of agreement of 0.03. The distribution of anti-Xa measurements in patients with and without clinically relevant drug levels (30, 50, and 100  $\mu$ gL<sup>-1</sup>) is shown in **Figure 3**.

Sensitivity analyses considering samples below  $3.34 \text{ U ml}^{-1}$  only yielded the following results similar to the values mentioned above:  $r_s$  overall 0.97 (95%CI 0.96, 0.97),  $r_s$  rivaroxaban 0.97 (0.96, 0.97),  $r_s$  apixaban 0.96 (0.95, 0.97), and  $r_s$  edoxaban 0.96 (0.93, 0.98). The Deming regression slope was 0.01 (overall; 95% CI 0.01, 0.01), 0.01 (rivaroxabar; 0.01, 0.01), 0.01 (apixabar; 0.01, 0.01), and 0.01 (edoxabar; 0.01, 0.01). The Y-intercept was -0.01 (overall; -0.05, 0.03), -0.03 (rivaroxabar; -0.08, 0.01), -0.01 (apixabar; -0.06, 0.04), and 0.04 (edoxabar; -0.05, 0.13).

# Diagnostic Accuracy Regarding Clinically Significant Drug Levels

In the derivation dataset (n = 422), the area under the ROC curve was 0.99 for all clinically relevant drug concentrations (**Figure 4**, panel A); 95% CI were 0.977 to 0.997 in case of 30 µg L<sup>-1</sup>, 0.981–0.995 in 50 µg L<sup>-1</sup>, and 0.978–0.994 in 100 µg L<sup>-1</sup>. Drug levels above 30 µg L<sup>-1</sup> were detected with a sensitivity of 96.9% (95% CI 94.3–98.5; cut-off value 0.26 U mL<sup>-1</sup>), drug levels above 50 µg L<sup>-1</sup> with a sensitivity of 95.2% (95% CI 92.0–97.4; cut-off value 0.40), and drug levels above 100 µg L<sup>-1</sup> with a sensitivity of 96.2% (95% CI 92.4–98.5; cut-off value 0.77). Specificities were 92.2% (85.1–96.6), 90.6% (84.7–94.8), 88.1% (83.3–92.0), respectively. Confusion matrices are given in **Figure 4**, panel B. In the sensitivity analyses, these measures did not differ significantly.

## **Internal Validation**

In the verification dataset (n = 423), the area under the ROC curve was 0.98 for both the cut-offs 30 (0.971–0.990) and 50 (0.976–0.993) µg L<sup>-1</sup> and 0.99 (0.988–0.998) for 100 µg L<sup>1</sup> (**Figure 5**, panel A). Cut-off levels obtained in the derivation dataset were used to calculate sensitivities and specificities in the verification dataset (n = 423) as an internal validation (**Figure 5**, panel B). Drug levels above 30 µg L<sup>-1</sup> were detected with a sensitivity of 94.2% (95% CI 90.8–96.6; cut-off value 0.26 U mL<sup>-1</sup>), drug levels above 50 µg L<sup>-1</sup>



(85.1, 96.6), 90.6% (84.7, 94.8), 88.1% (83.3, 92.0), respectively.

with a sensitivity of 95.8% (95% CI 92.4–98.0; cut-off value 0.40), and drug levels above 100  $\mu$ g L<sup>-1</sup> with a sensitivity of 98.7% (95% CI 95.5–99.9; cut-off value 0.77). Specificities were 86.3% (79.2–91.7), 89.8% (84.5–93.7), 88.7% (84.2–92.2), respectively. In the sensitivity analyses, these measures did not differ significantly.

# DISCUSSION

We conducted a large multicenter cross-sectional study including 845 patients treated with rivaroxaban, apixaban, or edoxaban in clinical practice. Even though minor differences between drugs and some variability in the higher concentrations exist, the association between heparin-calibrated anti-Xa measurements and LC-MS/MS results was strong for all drugs. As an internal validation, clinically relevant drug levels were predicted correctly in the verification dataset using cut-off levels derived in the derivation dataset.

Even though this is the first study analyzing a UFHcalibrated anti-Xa assay in a large cohort, our results are essentially in-line with previous publications. Recently, we analyzed the accuracy of a LMWH-calibrated anti-Xa assay in the same cohort (n = 867) (21). Similarly, the accuracy was high and clinically relevant drug levels were predicted



correctly. Studt et al. (17) studied the accuracy and consistency of anti-Xa assays for rivaroxaban plasma concentration in 20 healthy individuals and found a high agreement with drug concentrations measured by LC-MS/MS. Another crosssectional study (25) showed a strong correlation between heparin anti-factor Xa activity and LC-MS/MS in 30 patients treated with rivaroxaban. In another retrospective study,

<b>TABLE 3</b>   Implementation of a single, heparin-calibrated anti-Xa assay for the   measurement of rivaroxaban, apixaban, and edoxaban drug concentrations.				
	Cut-off value	Regression equation		
Clinical threshold				
30 µg/L DOAC	0.26 U/mL			
50 μg/L DOAC	0.40 U/mL			
100 µg/L DOAC	0.77 U/mL			
Drug concentrations				
Rivaroxaban		129 × [U/mL] – 37		
Apixaban		124 × [U/mL] – 17		
Edoxaban		166 × [U/mL] – 16		

The cut-off values could be given, which would allow direct clinical decisions to be made. Alternatively, the regression equations can be used to calculate the plasma concentrations of each compound.

a high degree of correlation between a heparin-calibrated anti-Xa assay and LC-MS/MS was observed in 24 patients taking rivaroxaban or apixaban (26). Similar results were also observed in studies assessing the correlation of anti-Xa measurements with LC-MS/MS in spiked samples (27– 30). A high correlation between UFH-calibrated anti-Xa measurements and rivaroxaban-/apixaban-calibrated anti-Xa results was found in a study assessing 241 left-over samples (20). Besides, van Pelt and colleagues proposed a universal anti-Xa assay reporting the inhibitory effect rather than drug concentrations (18).

This study is associated with certain strengths and limitations. The sample size is much larger than previous studies (n = 845), which increases statistical power and precision. Thus, we were able to split the dataset into a derivation and verification dataset, facilitating internal validation of the results. Another strength is that our study was designed as a multicenter study (nine laboratories) conducted in clinical practice, thus constituting a representative population. Therefore, the results can be translated straightforwardly to clinical practice. Also, LC-MS/MS was used as a reference standard which is considered the most accurate technique to measure DOAC plasma levels (21, 31). As a limitation, only a single heparin-calibrated assay was studied and other heparin-calibrated anti-Xa assays may perform differently. Several previous studies using various study designs found considerable differences among reagents (19, 32, 33). Additionally, the number of edoxaban-treated patients was lower compared to rivaroxaban and apixaban-treated patients. However, we believe that this is compensated by the large sample size.

The question now arises as to how this test can be applied in daily practice. The cut-off values mentioned in Table 3 can be given so that clinical decisions can be prompted directly (Table 3). Alternatively, the regression equations can be used to calculate the plasma concentrations of each compound (Table 3). However, the second approach must be used with caution. Minor differences between drugs exist and a certain degree of variability can be observed in higher concentrations. Therefore, drug levels in high concentrations ( $\geq 150 \ \mu g L^{-1}$ ) can only be measured with limited accuracy and precision (34). Yet, several arguments can be raised in favor of the heparin-calibrated assay: (1) this variation is also present in drug-specific anti-Xa assays (17), (2) typical requirements in terms of correlation coefficients and ROC AUC are fulfilled, and (3) all clinically relevant cut-off thresholds are met with high accuracy. Current guidelines do not distinguish between normal-high and very-high drug levels concerning patient management.

Our results confirm that a single, heparin-calibrated anti-Xa assay accurately measures rivaroxaban, apixaban, and edoxaban drug concentrations and correctly predicts clinically relevant drug levels. Since heparin-calibrated anti-Xa assays are already available in many laboratories, determination of DOAC drug levels can be provided easily. This may foster the widespread implementation of anti-Xa assays to measure DOAC, thus improving care in patients taking these drugs. Future studies should confirm these results in other settings and using other heparin-calibrated assays.

## Conclusion

We report results of a large prospective study including patients treated with rivaroxaban, apixaban, or edoxaban in clinical practice. The association of heparin-calibrated anti-Xa measurements with DOAC drug concentrations was strong, and clinically relevant drug levels were predicted correctly. Our results represent a strong argument in favor of the potential application of a universal, heparin-calibrated anti-Xa assay to measure DOAC in clinical practice.

# DATA AVAILABILITY STATEMENT

The raw data supporting the conclusions of this article will be made available by the authors, without undue reservation.

# ETHICS STATEMENT

The studies involving human participants were reviewed and approved by the Kantonale Ethikkommission Bern. The patients/participants provided their written informed consent to participate in this study.

# **AUTHOR CONTRIBUTIONS**

TM analyzed the data, interpreted the results, and wrote the manuscript. J-DS, AM, LAI, PF, WW, AS, LG, BG, CB, LAs, UA, and TS collected data, and contributed to study design, protocol, and preparation of the manuscript. CB and UA contributed essential tools and reagents. MN designed the study, wrote the protocol, collected data, analyzed, and interpreted the data, and wrote the manuscript. All authors contributed to the article and approved the submitted version.

# FUNDING

This study was supported by a research grant of the Research Fund Haematology Cantonal Hospital Lucerne. MN was supported by a research grant of the Swiss National Science Foundation (#179334). Implementation of the LC-MS/MS measurements was supported by the Gottfried & Julia Bangerter-Rhyner-Stiftung. This study was supported by a research grant of Siemens Healthineers. These funders had no role in study design, data collection and analysis, the decision to publish, or manuscript preparation. We thank the following companies for the provision of reagents and/or pure substances: Bayer Healthcare AG, Bristol-Myers Squibb, and Daiichi Sankyo. These companies had no role in study design, data collection and analysis, the decision to publish, or manuscript preparation.

## ACKNOWLEDGMENTS

We thank Evelyne Giabbani for organizing the sample collection, Marc Gisi, and Vincent Benites for doing the anti-Xa measurements, Gabriela Mäder for doing the LC-MS/MS measurements, Sarah Phillip for technical support, and Martin Fiedler for providing the infrastructure at Inselspital.

## REFERENCES

- Barnes GD, Lucas E, Alexander GC, Goldberger ZD. National trends in ambulatory oral anticoagulant use. Am J Med. (2015) 128:1300–5.e2. doi: 10.1016/j.amjmed.2015.05.044
- van den Heuvel JM, Hovels AM, Buller HR, Mantel-Teeuwisse AK, de Boer A, Maitland-van der Zee AH. NOACs replace VKA as preferred oral anticoagulant among new patients: a drug utilization study in 560 pharmacies in the Netherlands. *Thromb J.* (2018) 16:7. doi: 10.1186/s12959-017-0156-y
- Sauter TC, Eberle B, Wuillemin WA, Thiele T, Angelillo-Scherrer A, Exadaktylos AK, et al. How I manage patients with anticoagulation-associated bleeding or urgent surgery. *Swiss Med Wkly.* (2018) 148:w14598. doi: 10.4414/ smw.2018.14598
- Cuker A. Laboratory measurement of the non-vitamin K antagonist oral anticoagulants: selecting the optimal assay based on drug, assay availability, and clinical indication. *J Thromb Thrombolysis*. (2016) 41:241–7. doi: 10.1007/ s11239-015-1282-7
- Cuker A, Siegal D. Monitoring and reversal of direct oral anticoagulants. *Hematology Am Soc Hematol Educ Program.* (2015) 2015:117–24. doi: 10.1182/ asheducation.v2015.1.117.3916182
- Kaserer A, Kiavialaitis GE, Braun J, Schedler A, Stein P, Rossler J, et al. Impact of rivaroxaban plasma concentration on perioperative red blood cell loss. *Transfusion*. (2020) 60:197–205. doi: 10.1111/trf.15560
- Spahn DR, Bouillon B, Cerny V, Duranteau J, Filipescu D, Hunt BJ, et al. The European guideline on management of major bleeding and coagulopathy following trauma: fifth edition. *Crit Care.* (2019) 23:98. doi: 10.1186/s13054-019-2347-3
- Levy JH, Ageno W, Chan NC, Crowther M, Verhamme P, Weitz JI, et al. When and how to use antidotes for the reversal of direct oral anticoagulants: guidance from the SSC of the ISTH. J Thromb Haemost. (2016) 14:623–7. doi: 10.1111/jth.13227
- Lindhoff-Last E. Direct oral anticoagulants (DOAC) management of emergency situations. *Hamostaseologie*. (2017) 37:257–66. doi: 10.5482/ HAMO-16-11-0043
- Douxfils J, Ageno W, Samama CM, Lessire S, Ten Cate H, Verhamme P, et al. Laboratory testing in patients treated with direct oral anticoagulants: a practical guide for clinicians. *J Thromb Haemost.* (2018) 16:209–19. doi: 10.1111/jth.13912
- Moner-Banet T, Alberio L, Bart PA. Does one dose really fit all? On the monitoring of direct oral anticoagulants: a review of the literature. *Hamostaseologie*. (2020) 40:184–200. doi: 10.1055/a-1113-0655
- Kaserer A, Schedler A, Jetter A, Seifert B, Spahn DR, Stein P, et al. Risk factors for higher-than-expected residual rivaroxaban plasma concentrations in real-life patients. *Thromb Haemost.* (2018) 11:808–17. doi: 10.1055/s-0038-1639585
- Xu Y, Schulman S, Dowlatshahi D, Holbrook AM, Simpson CS, Shepherd LE, et al. Healthcare resource utilization and costs among patients with direct oral anticoagulant or warfarin-related major bleeding. *Thromb Res.* (2019) 182:12–9. doi: 10.1016/j.thromres.2019.07.026
- Ebner M, Peter A, Spencer C, Hartig F, Birschmann I, Kuhn J, et al. Point-ofcare testing of coagulation in patients treated with non-vitamin K antagonist oral anticoagulants. *Stroke*. (2015) 46:2741–7. doi: 10.1161/STROKEAHA.115. 010148
- 15. Hawes EM, Deal AM, Funk-Adcock D, Gosselin R, Jeanneret C, Cook AM, et al. Performance of coagulation tests in patients on therapeutic doses of dabigatran: a cross-sectional pharmacodynamic study based on peak and

## SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at: https://www.frontiersin.org/articles/10.3389/fcvm. 2022.817826/full#supplementary-material

**Supplementary Figure 1** | A modified Bland–Altman plot using ratios to determine a possible systematic bias over the range of measurements. The bias and the upper and lower limit of agreement is shown.

trough plasma levels. J Thromb Haemost. (2013) 11:1493–502. doi: 10.1111/ jth.12308

- Samama MM, Meddahi S, Samama CM. Pharmacology and laboratory testing of the oral Xa inhibitors. *Clin Lab Med.* (2014) 34:503–17. doi: 10.1016/j.cll. 2014.06.009
- Studt JD, Alberio L, Angelillo-Scherrer A, Asmis LM, Fontana P, Korte W, et al. Accuracy and consistency of anti-Xa activity measurement for determination of rivaroxaban plasma levels. *J Thromb Haemost.* (2017) 15:1576–83. doi: 10.1111/jth.13747
- van Pelt LJ, Lukens MV, Testa S, Chatelain B, Douxfils J, Mullier F. The DaXainhibition assay: a concept for a readily available, universal aXa assay that measures the direct inhibitory effect of all anti-Xa drugs. *Thromb Res.* (2018) 168:63–6. doi: 10.1016/j.thromres.2018.04.024
- Sabor L, Raphael M, Dogne JM, Mullier F, Douxfils J. Heparin-calibrated chromogenic anti-Xa assays are not suitable to assess the presence of significant direct factor Xa inhibitors levels. *Thromb Res.* (2017) 156:36–8. doi: 10.1016/j.thromres.2017.05.024
- Maier CL, Asbury WH, Duncan A, Robbins A, Ingle A, Webb A, et al. Using an old test for new tricks: measuring direct oral anti-Xa drug levels by conventional heparin-calibrated anti-Xa assay. *Am J Hematol.* (2019) 94:E132-4. doi: 10.1002/ajh.25434
- 21. Willekens G, Studt JD, Mendez A, Alberio L, Fontana P, Wuillemin WA, et al. A universal anti-Xa assay for rivaroxaban, apixaban, and edoxaban measurements: method validation, diagnostic accuracy, and external validation. *Br J Haematol.* (2021) 193:1203–12. doi: 10.1111/bjh.17470
- 22. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline. CLSI Document H21-A5. 5th ed. Wayne, PA: Clinical and Laboratory Standards Institute (2008).
- Horber S, Lehmann R, Peter A. Evaluation of the Atellica COAG 360 coagulation analyzer in a central laboratory of a maximum care hospital. *Int J Lab Hematol.* (2020) 42:28–36. doi: 10.1111/ijlh.13119
- Pollock MA, Jefferson SG, Kane JW, Lomax K, MacKinnon G, Winnard CB. Method comparison–a different approach. Ann Clin Biochem. (1992) 29(Pt 5):556–60. doi: 10.1177/000456329202900512
- Gosselin RC, Francart SJ, Hawes EM, Moll S, Dager WE, Adcock DM. Heparin-calibrated chromogenic anti-Xa activity measurements in patients receiving rivaroxaban: can this test be used to quantify drug level? *Ann Pharmacother.* (2015) 49:777–83. doi: 10.1177/1060028015578451
- Beyer J, Trujillo T, Fisher S, Ko A, Lind SE, Kiser TH. Evaluation of a heparincalibrated antifactor Xa assay for measuring the anticoagulant effect of oral direct Xa inhibitors. *Clin Appl Thromb Hemost.* (2016) 22:423–8. doi: 10.1177/ 1076029616629759
- Asmis LM, Alberio L, Angelillo-Scherrer A, Korte W, Mendez A, Reber G, et al. Rivaroxaban: quantification by anti-FXa assay and influence on coagulation tests: a study in 9 Swiss laboratories. *Thromb Res.* (2012) 129:492–8. doi: 10.1016/j.thromres.2011.06.031
- Samama MM. Which test to use to measure the anticoagulant effect of rivaroxaban: the anti-factor Xa assay. J Thromb Haemost. (2013) 11:579–80. doi: 10.1111/jth.12165
- Bardy G, Fischer F, Appert A, Baldin B, Steve M, Spreux A, et al. Is antifactor Xa chromogenic assay for rivaroxaban appropriate in clinical practice? Advantages and comparative drawbacks. *Thromb Res.* (2015) 136:396–401. doi: 10.1016/j.thromres.2015.05.015
- Schmitz EM, Boonen K, van den Heuvel DJ, van Dongen JL, Schellings MW, Emmen JM, et al. Determination of dabigatran, rivaroxaban and apixaban

by ultra-performance liquid chromatography – tandem mass spectrometry (UPLC-MS/MS) and coagulation assays for therapy monitoring of novel direct oral anticoagulants. *J Thromb Haemost.* (2014) 12:1636–46. doi: 10.1111/jth.12 702

- 31. Schellings MW, Boonen K, Schmitz EM, Jonkers F, van den Heuvel DJ, Besselaar A, et al. Determination of dabigatran and rivaroxaban by ultra-performance liquid chromatography-tandem mass spectrometry and coagulation assays after major orthopaedic surgery. *Thromb Res.* (2016) 139:128–34. doi: 10.1016/j.thromres.2016.01.012
- Smahi M, De Pooter N, Hollestelle MJ, Toulon P. Monitoring unfractionated heparin therapy: lack of standardization of anti-Xa activity reagents. *J Thromb Haemost.* (2020) 18:2613–21. doi: 10.1111/jth.14969
- Hollestelle MJ, van der Meer FJM, Meijer P. Quality performance for indirect Xa inhibitor monitoring in patients using international external quality data. *Clin Chem Lab Med.* (2020) 58:1921–30. doi: 10.1515/cclm-2020-0130
- 34. Margetic S, Celap I, Delic Brkljacic D, Pavlovic N, Supraha Goreta S, Kobasic I, et al. Chromogenic anti-FXa assay calibrated with low molecular weight heparin in patients treated with rivaroxaban and apixaban: possibilities and limitations. *Biochem Med (Zagreb)*. (2020) 30:010702. doi: 10.11613/BM.2020. 010702

**Conflict of Interest:** MN reports research grants from Bayer Healthcare, outside of the submitted work, lecture honoraria from Bayer Healthcare, and Daiichi Sankyo. LAI reports research grants from Bayer, CSL-Behring, Novartis, Novo Nordisk, Roche, Sobi, and Takeda. WW reports research grants from Bayer Healthcare, BMS-Pfizer, Daiichi Sankyo and Sanofi, and honoraria for participating in scientific advisory boards from Bayer, Pfizer, and from Alexion Pharma GmbH, all outside the submitted work. J-DS reports lecture fees and

advisory honoraria from Bayer Healthcare, Pfizer, Takeda, Siemens, and Sanofi. TS holds an endowed professorship supported by the Touring Club Switzerland. BG reports non-financial support and funding for accredited continuing medical education program from Axonlab, and from Thermo Fisher Scientific, during the conduct of the study; personal fees and funding for accredited continuing medical education program from Alnylam, grants, personal fees and funding for accredited continuing medical education program from Pfizer, funding for accredited continuing medical education program from Bayer, Bristol Myers Squibb, Daiichi Sankyo, Takeda, Octapharma, SOBI, Janssen, Novo Nordisk, Mitsubishi Tanabe Pharma, outside the submitted work.

The remaining authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

**Publisher's Note:** All claims expressed in this article are solely those of the authors and do not necessarily represent those of their affiliated organizations, or those of the publisher, the editors and the reviewers. Any product that may be evaluated in this article, or claim that may be made by its manufacturer, is not guaranteed or endorsed by the publisher.

Copyright © 2022 Meihandoest, Studt, Mendez, Alberio, Fontana, Wuillemin, Schmidt, Graf, Gerber, Amstutz, Bovet, Sauter, Asmis and Nagler. This is an openaccess article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.